



## INNOVATOR X<sup>®</sup> POST-OP ELBOW

Benefiting from our proven, patented technologies, Innovator X, with extending struts, can help improve patient outcomes. With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, Innovator X was designed to increase patient satisfaction and compliance.

- Patented Form Fit<sup>®</sup> 3-dimensional padding and Lycra<sup>®</sup> lining improve fit and comfort
- Comfort pad at the neck and shoulder
- Optional arm bar accessory kit provides anti-supination/pronation control during use



- 1 OneTouch button enables easy release for slide-to-size length adjustments
- 2 Easy-to-use Innovator hinge makes flexion and extension adjustments simple
- 3 Sensil<sup>®</sup> silicone padding on straps reduces migration for an improved fit and increased stability
- 4 Malleable aluminum cuffs allow customization for increased fit and support

## INDICATIONS

- Post elbow surgery: elbow reconstruction, ligament and tendon repairs, surgical repair of a fracture and biceps repair
- Stable fractures
- Mild to severe soft tissue injuries

PART#	DESCRIPTION
W-823101	Left
W-823100	Right
W-A82310	Arm bar kit

Patent #: 5,334,135; 7,037,287  
Other U.S. and foreign patents pending

Simply press, rotate and release flexion and extension setting